 **TALASSI**

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**SERIE
FILTERBOX**

FILTRES
A CHARBON ACTIF

**SERIE
FILTERBOX**

ACTIVATED
CARBON FILTERS

**SERIES
FILTERBOX**

AKTIVKOHLEFILTER

**SERIE
FILTERBOX**

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PVC ACTIVATED CARBON FILTERS FILTERBOX

The control of atmospheric emissions involves both the comfort and the ecological aspect of working environments.

Respect for the environment, the EEC directives, and the well-being of the people who work in laboratories are all factors that require the use of systems which hold back acids or solvents present in chemical processing discharges, with a resulting air purification.

The **FilterBox** PVC activated carbon filters are the answer provided by **Plastifer** to the problem of air filtering.

Main characteristics:

Great adsorption capacity, due to the activated carbon plates especially designed to filter fumes in laboratories.

Optimum aeraulic efficiency, due to the special pipe fittings that optimise air distribution inside the filter, reducing pressure losses.

Completely manufactured in acid-proof PVC. All its components are resistant to chemical agents.

High-performance pre-filter

It protects the activated carbon from dust, which increases performance and durability; it has an anticorrosive injection-moulded housing that makes it easy to be replaced.

It can be installed outdoors, thanks to the PVC cover that protects it against rainfall and its polypropylene screws.

Easy and rapid installation combines with reduced sizes.

Safe maintenance operations can be performed thanks to the compact carbon plates that has been designed to avoid the risk of breathing in noxious powders during replacement operations.

Cost-effectiveness: the separation of the pre-filter from the carbon plates enables the pre-filter to be cleaned without having to replace the carbon plates as well.

Plastifer recommends installing **FilterBox** combined with the series of **Polyfan** anti-acid aspiration units.

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CARBON PLATES PERFORMANCE

Based on a strong partnership with the world's leading activated carbon producer, we have tested and selected more than 150 types of activated carbon to find the one that most appropriately adsorbs fumes in laboratories.

The highly efficacious Plastifer plates are made with highly efficient activated carbons (the same ones used in industrial gas masks) due to the distribution of pores that guarantees perfect adsorption of the compounds.

The Regular Plastifer plate is made with activated carbons approved for class A and AX.

The RBAA Plastifer plate is made with activated carbons approved for class E.

The Regular and RBAA Plastifer activated carbon plates are products exclusively manufactured for Plastifer.

WHAT IS ACTIVATED CARBON:

Activated carbon is a microporous granular carbon that can be produced with different raw material such as peat, wood or fossil coal.

The activation process, carried out with vapour or chemical agents, produces millions of pores of molecular dimension that make up a great surface area with a ratio probably larger than 2,000 m² per gram of activated carbon.

Due to the Van Der Waals force, the carbon atoms present on the internal surface of the activated carbon attract the neighbouring gas molecules.

The intensity of such phenomenon is related to the molecular structure of the suctioned substances; some molecules are strongly attracted whereas others are not. Due to this mechanism, activated carbon can be used to eliminate one or more components of a mixture.

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CARBON PLATES PERFORMANCE

DATA TO CONSIDER WHEN CHOOSING A FILTER:

- Airflow to be treated.
- Composition and concentration of the air to be treated.
- Required performance.
- Maximum temperature 60°.
- Relative humidity lower than 70 %.

CAUSES THAT REDUCE CARBON PERFORMANCE:

Temperature reduces the carbon adsorption capacity. The lower the temperature, the higher the carbon adsorption capacity.

Humidity has a negative influence on the carbon adsorption capacity since it can withstand a maximum relative humidity equal to 70%.

Dust is another enemy of carbon because it blocks its pores and the spaces between the different cylinders consequently reducing performance and increasing pressure losses.

USEFUL RECOMMENDATION:

The **contact time** between the air and the activated carbon is an important factor to consider; on the basis of the required result we suggest values that range between 0,05 sec. and 0,2 sec.

The activated carbon **compactness** is a fundamental aspect that must be observed in order to achieve effective results.

Plates and panels that are not compacted or that have great granulometric values cannot render good results.

The **installation** of the activated carbon filter must enable it to work in depression.

The use of very small activated carbon filters should be discouraged.

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